

⚠ Safety Instructions

WARNING: Not to be used as a visual public mode alarm notification

1. Read and understand all instructions.
2. Follow all warnings and instructions on the Simplicity™ Telephone Ring Signaler LTD.
3. For indoor use only.
4. Unplug the Signaler from the wall outlet, telephone and lamp before cleaning. Do not use liquid or aerosol cleaners. Use a damp cloth for cleaning.
5. Do not use the Signaler near water, for example, near a bathtub, in a wet basement, or near a swimming pool.
6. The Signaler should be operated only from the type of power source indicated on the marking label.
7. Do not allow anything to rest on the cords of the Signaler or of the lamp. Do not locate the Signaler where people can

FCC Information

The Telephone Ring Signaler LTD has been tested and found to comply with the specifications for a Class B digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates and uses radio frequency energy and it not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the Signaler and the receiver.
- Connect the Signaler to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA. A label is located on the underside of the Telephone Ring Signaler LTD containing, among other information, an FCC registration number, DAKUSA-35957-OT-N. If requested, this number must be provided to the telephone company. A plug and jack used to connect this equipment to the premises wiring and telephone network must comply with the applicable FCC Part 68 rules and requirements adopted by the ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. See installation instructions for details.

The RBN is used to determine the quantity of devices which may be connected to the telephone line. Excessive RBNs on the line may result in the devices not ringing in response to an incoming call. In most, but not all, areas the sum of RBNs should not exceed five (5, 0). To be certain

8. Do not overload wall outlets and extension cords. Walk on the cords.
9. To reduce the risk of electrical shock, do not disassemble the Signaler. Take it to a qualified service person when service or repair work is required. Opening the unit or incorrect re-assembly can expose you to dangerous voltages or other risks. The Signaler does not have any user-serviceable parts.
10. Unplug the Signaler from the wall outlet and refer servicing to qualified service personnel under the following conditions: (a) when the plug is damaged; (b) if liquid has been spilled into the Signaler; (c) if the Signaler has been exposed to rain or water; (d) if the Signaler does not operate normally; (e) if the Signaler has been dropped and/or damaged; (f) if the Signaler exhibits a distinct change in performance.

of the number of devices that may be connected to a line, as determined by the total RBN, contact the local telephone company. The RBN for this product is 0,7B.

Should you experience trouble with this telephone equipment, please contact: Ultratec National Service Center, 5901 Research Park Blvd., Madison, WI 53719, email: service@ultratec.com.

For repair or warranty information, please contact Ultratec Customer Service at (608) 238-5400 (V/TTY). If the equipment is causing harm to the telephone network, the telephone company may request that you disconnect the equipment until the problem is resolved.

This equipment cannot be used on public coin phone service provided by the telephone company. Connection to party line service is subject to state tariffs.

If the Signaler causes harm to the telephone network, the telephone company will notify you in advance that temporary discontinuance of service may be required. But, if advance notice is not practical, the telephone company will notify you as soon as possible. Also, you will be advised of your right to file a complaint with the FCC if you believe it is necessary.

Occasionally, your telephone company may make changes in its facilities, operation, or procedures that could affect the operation of your equipment. If so, you will be given advance notice of the change to give you an opportunity to maintain uninterrupted service.

If your home has specially wired alarm equipment connected to the telephone line, ensure the installation of the Signaler does not disable it. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.

Other than replacing the batteries, the Telephone Ring Signaler LTD does not have any user-serviceable parts. Modification or changes to the Signaler not expressly approved by Ultratec, Inc. can void your authority to operate the Signaler.


11. Adjust only those controls that are covered by the operating instructions. Improper adjustment of other controls may result in damage.

Where to Write

Your comments can help improve our products, services and instruction manuals. If you have an idea or suggestion, or if you need more information about your Simplicity units, contact:

Ultratec, Inc.
450 Science Dr.
Madison, WI 53711
(800) 482-2424 (Voice/TTY)
Email: service@ultratec.com
www.ultratec.com

Using your *Simplicity*™ Telephone Ring Signaler LTD



©2001, 2002 Ultratec, Inc.
 July 2002



305-012201

Two-Year Limited Warranty

Ultratec, Inc. offers a two-year warranty for the *Simplicity*™ Telephone Ring Signaler LTD. Warranty service is available through Ultratec. If your product needs service, return it to:

Ultratec National Service Ctr.
5901 Research Park Blvd.
Madison, WI 53719

Warranty Terms

Ultratec shall repair or replace a defective product, at Ultratec's option, free of charge if returned to Ultratec, freight prepaid, within two years after the product has been shipped to Purchaser, provided, however, that Purchaser shall give Ultratec written notice within said time period of the fact of such defective product and shall also indicate in writing the nature of the defect upon returning the defective product to Ultratec.

This warranty covers defects in materials and workmanship only. It does not cover damage caused by accidents, defects, or damage resulting from misuse, abuse, or disregard for reasonable care, or damage caused by subjecting the product to unspecified electrical voltage or attaching inappropriate devices (including, but not limited to, damage caused by lightning and connecting the product to digital telephone systems).


Other than as stated herein, there are no other warranties, whether by sample or model, by description or quantity, whether expressed or implied, and there are no implied warranties of fitness for a particular purpose or of merchantability.

This warranty shall be limited to the foregoing obligations of Ultratec and specifically excludes any incidental or consequential damages resulting from defects in the product.

In the case that the product specified above is shipped to an Ultratec service center for warranty service, the owner is responsible for payment of shipping charges. Because damage during shipment is not covered by this warranty, the product should be insured for shipment by the owner.

Setup

1. Plug the lamp into the Signaler plug.
 2. Plug the telephone line cord into any unused telephone outlet and plug the other end into the Signaler. To plug the telephone into the same outlet, use a "Y-jack" (also called a duplex jack).
 3. Plug the Signaler into an electrical outlet.
 4. Using the switch on the lamp, turn the lamp on. The switch on the lamp must always remain on to signal sounds.
- If you want to turn the lamp on and off for other purposes such as reading, use the Signaler **On/Off Switch**.




- ## Troubleshooting
- ### Transmitting Problems
- If the signal lamp doesn't flash when the phone rings, check the following:
- If the red **Transmit Light** on the Signaler does not blink when the phone rings, the Signaler is not properly connected.
 - Be sure the Signaler is properly plugged into a working electrical outlet and a working telephone line.
 - Be sure that, if you electrical outlet is controlled by a wall switch, the switch is on.
 - Do not plug the Signaler into a power strip.
 - If there is more than one *Simplicity* System operating in your building, they may interfere with each other.
- ### Lamp Problems
- If the red **Transmit Light** on the Signaler blinks when the phone rings, the signaler is properly connected and the problem is in the lamp
- connection. Check the following:
- Be sure the lamp is plugged into the Signaler.
 - Be sure the lamp switch is on, and the lightbulb is working.

Simplicity™ System


Other units that work with the Telephone Ring Signaler LTD:



Simplicity™ Remote Receiver LR



Simplicity™ Remote Receiver LRD



Simplicity™ Remote Receiver SRD

Specifications

Physical Dimensions

Size: 3.75 x 2.75 x 1.25 in.
Weight: 11 oz.

Connections

RJ11 Outlet
120 Volt AC Outlet
300 Watts maximum load

Controls

Lamp On/Off Switch
Transmit Light

Approvals

FCC approved

FCC Reg. No.: D8K USA-35957-OT-N
Safety: CSA NRTL/C both US and Canada

Canada

Industry Canada: 721 11286A

RBN 0.7B USA 0.5 Canada